


# Module: Shimmer GSR

## Base & Advanced Packages

Shimmer3 GSR Wireless Device

Pre-requisite: iMotions Core License



 **IMOTIONS®**

[sales@imotions.com](mailto:sales@imotions.com)



# Shimmer GSR Module

Enables iMotions software to connect to the Shimmer GSR device. The Module is divided in Basic & Advanced Packages



## Shimmer3 GSR wireless device:

Shimmer3 GSR device provides real-time physiological data collection and wireless streaming with high-quality, scientifically reliable data. The Shimmer3 monitors & wirelessly transmit skin conductivity GSR (Galvanic Skin Response) and optical heart rate between two reusable electrodes attached to two fingers of one hand and heart rate with a third electrode. Measured as increases in perspiration at the surface of the skin, a byproduct of sympathetic nervous system activity, GSR represents a measure of arousal or interest.

FEATURES	Basic Module	Advanced Module
GSR Signals	✓	✓
Optical Heart Rate Stream	✓	✓
Accelerometer Signals	✓	✓
Stimuli Presentation	✓	✓
Real Time Synchronization	✓	✓
Raw Data Export	✓	✓
Visualizations Export	✓	✓
Live / Post Markers	✓	✓
Individual & Aggregated View	✓	✓
Signal Quality Assurance	✓	✓
Battery Level Report	✓	✓
GSR Peak Aggregation		✓
Multi-Shimmer Integration		✓

Some users of iMotions Shimmer GSR Module




The Shimmer GSR Module has the iMotions Core License as pre-requisite, which allows you to:



### **Present all Type of Stimuli**

Present images, videos, websites, games, screen & scene recordings, real world products and surveys.



### **Create Sophisticated Studies**

Full flexibility to design any study setup with randomizations, block designs, test plans, group rotations and more.



### **Real-Time Synchronization**

GSR data, stimuli and any other sensor data streams are real time synchronized.




### **Export Raw Data**

Export all collected GSR data in sync with stimuli and other sensors in .txt format for easily import into MatLab and other statistical programs.



### **Get Quality Assurance**

Monitor data collection quality at any given time to ensure the validity of the studies.



### **Create Live or Post Markers**

Mark important happenings during data collection or in replay mode to facilitate the analysis.



### **View Individual & Aggregated**

Get visualisations whether individually or aggregated in any type of segmentation and export them as images and/or videos




### **Execute Studies in a Lab and/or Mobile Environments**

iMotions Shimmer GSR Module is a wireless system that allows for any study design, whether this is a controlled environment setup or a mobile scenario where respondents move freely in a space.



# The Shimmer GSR module delivers visualizations of collected GSR data in sync with any type of stimuli



Synchronized GSR raw data is also available by the click of a button

The screenshot shows an Excel spreadsheet with the following columns:

- StimulusName
- SlideType
- Timestamp
- MediaTime
- Cal Timestamp (SEE8)
- Cal VSenseBatt (SEE8)
- Cal Heart Rate [BPM] (SEE8)
- Cal GSR (SEE8)
- GSR Valid Value (SEE8)

StimulusName	SlideType	Timestamp	MediaTime	Cal Timestamp (SEE8)	Cal VSenseBatt (SEE8)	Cal Heart Rate [BPM] (SEE8)	Cal GSR (SEE8)	GSR Valid Value (SEE8)
BASLINE	TestImage	20151030_104845058	0	203769.5313	3831.501832	146	15.8656	VALID
BASLINE	TestImage	20151030_104845063	0	203779.2969	3831.501832	146	15.8605	VALID
BASLINE	TestImage	20151030_104845063	0	203789.0625	3834.432234	146	15.8656	VALID
BASLINE	TestImage	20151030_104845073	0	203798.8281	3834.432234	146	15.8605	VALID
BASLINE	TestImage	20151030_104845083	0	203808.5938	3815.384615	146	15.8656	VALID
BASLINE	TestImage	20151030_104845093	0	203818.3594	3832.967033	146	15.8605	VALID
BASLINE	TestImage	20151030_104845103	0	203828.125	3830.03663	146	15.8605	VALID
BASLINE	TestImage	20151030_104845119	0	203837.8906	3828.571429	146	15.8554	VALID
BASLINE	TestImage	20151030_104845119	0	203847.6563	3822.710623	146	15.8554	VALID
BASLINE	TestImage	20151030_104845128	0	203857.4219	3809.52381	146	15.8554	VALID
BASLINE	TestImage	20151030_104845139	0	203867.1875	3830.03663	146	15.8554	VALID
BASLINE	TestImage	20151030_104845154	0	203876.9531	3832.967033	146	15.8503	VALID
BASLINE	TestImage	20151030_104845163	0	203886.7188	3824.175824	146	15.8503	VALID
BASLINE	TestImage	20151030_104845170	0	203896.4844	3830.03663	146	15.8503	VALID
BASLINE	TestImage	20151030_104845188	0	203906.25	3810.989011	146	15.8503	VALID
BASLINE	TestImage	20151030_104845189	0	203916.0156	3821.245421	146	15.8503	VALID
BASLINE	TestImage	20151030_104845200	0	203925.7813	3835.897436	146	15.8452	VALID
BASLINE	TestImage	20151030_104845227	0	203935.5469	3830.03663	146	15.8452	VALID
BASLINE	TestImage	20151030_104845228	0	203945.3125	3832.967033	146	15.8452	VALID
BASLINE	TestImage	20151030_104845232	0	203955.0781	3832.967033	146	15.8452	VALID
BASLINE	TestImage	20151030_104845243	0	203964.8438	3813.919414	146	15.8452	VALID
BASLINE	TestImage	20151030_104845244	0	203974.6094	3830.03663	146	15.8452	VALID
BASLINE	TestImage	20151030_104845254	0	203984.375	3821.245421	146	15.8401	VALID
BASLINE	TestImage	20151030_104845267	0	203994.1406	3831.501832	146	15.835	VALID
BASLINE	TestImage	20151030_104845277	0	204003.9063	3822.710623	146	15.835	VALID
BASLINE	TestImage	20151030_104845287	0	204013.6719	3812.454212	146	15.835	VALID
BASLINE	TestImage	20151030_104845297	0	204023.4375	3832.967033	146	15.8401	VALID
BASLINE	TestImage	20151030_104845309	0	204033.2031	3835.897436	146	15.835	VALID
BASLINE	TestImage	20151030_104845322	0	204042.9688	3822.710623	146	15.8299	VALID






# Shimmer GSR Advanced Package

Enables to automatically detect peaks in the GSR data of participants & aggregate them in any customized segmentation. It also allows the integration of up to 7 Shimmers whether they are GSR, EMG or ECG devices collecting data at the same time.

## 1. Automatic Peak Detection Aggregation



*Image shows visualizations of stimulus (video) and Peak Detection aggregation in time line. The numbers to the right indicate the number of individuals who displayed a peak response at a specific time.*


## 2. Multi-Shimmer Integration



This feature enables synchronization of a number of new data collection options such as GSR, EMG, and/or ECG signals, picked up simultaneously.

# Optionally combine GSR measurements with:

Eye Tracking, Facial Expression Recognition, EEG, ECG, EMG & Surveys



## Standard Lab view with combination of sensors





## Shimmer3 GSR

### Package contains:

- Hardshell box
- Bluetooth device
- Optical heart rate sensor (ear clip)
- 2 GSR sensors (cables + velcro straps + electrodes)
- Printed instruction / setup manual
- Hand strap

Technical Specifications	
Current Draw	60uA
Measurement Range	10k-4.7 MQ (.2uS - 100uS) +/- 10% 22k-680KQ (1.5-45uS) +/- 3%
Frequency Range	DC-15.9 Hz
Input Protection	RF/EMI filtering, Current limiting
Connections	Input 1, Input 2. All connections via Hospital-Grade 1.5 mm Touchproof IEC/EN 60601-1 DIN42-802 jacks
Dimensions	53 mm x 32 mm x 23 mm
Weight	32 g
Processing	MSP 430 microcontroller (8mHz, 16Bit)
Communication	Bluetooth - RN42, 802.15.4 radio - TI CC2420
Storage	Integrated 2GB microSD card slot
Battery	450mAh rechargeable Li-ion
Integrated 3 Axis Accel	Freescale MMA7361
Accel Range	+/- 1.5g-6g

## Relevant links:

Shimmer's technology has been referenced in over 100 peer review journal articles. Due to the technology's highly flexible nature, published papers include a broad range of application areas, from wireless sensor network research papers to falls detection, validating Shimmer's technology for a number of uses and against leading gold standard equipment.

1. [Event detection in an assisted living environment.](#)
2. [Real-time joint angle measurement using the Shimmer wireless sensor platform](#)
3. [SHIMMER™ – A Wireless Sensor Platform for Noninvasive Biomedical Research](#)
4. [SHIMMER and BioMOBIUS: A Health Research Platform](#)
5. [Automated Real-Time Atrial Fibrillation Detection on a Wearable Wireless Sensor Platform](#)
6. [A wireless framework for movement activity monitoring of sprinters](#)
7. [Research project: Shimmer Sensors for activity monitoring in post stroke upper limb rehabilitation](#)
8. [A WBAN for Human Movement Kinematics and ECG Measurements](#)
9. [Evaluation of Falls Risk in Community-Dwelling Older Adults Using Body-Worn Sensors](#)

For more publications visit: <https://delicious.com/ShimmerResearch>